


Anti-Ppan Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nq55

Item Number: CM0029719



Introduction

Rabbit polyclonal antibody targeting human PPAN (Suppressor of SWI4 1 homolog, Brix domain-containing protein 3). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse. Suitable for cell growth and nucleolar research.

[Learn More](#)

Parameter	Value
Product Name	PPAN Rabbit pAb
Remarks/Alias	SSF; SSF1; SSF2; BXDC3; SSF-1; PPAN
Species	Human
Gene ID (Human)	56342
Gene ID	56342
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PPAN (NP_064615.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQ55
Protein Weight	53kDa
Shipping	Ice bag