

Anti-Ppap2A Polyclonal Antibody For Wb, Elisa - O14494

Item Number: CM0030161



Introduction

PPAP2A Rabbit pAb targets human phospholipid phosphatase 1 (LPP1/PAP2a). Validated for WB and ELISA; reacts with human, mouse, rat. Ideal for lipid phosphate signaling and phosphatase research.

[Learn More](#)

Parameter	Value
Product Name	PPAP2A Rabbit pAb
Remarks/Alias	LPP1; PAP2; LLP1a; PAP-2a; PPAP2A
Species	Human
Gene ID (Human)	8611
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PPAP2A (NP_003702.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14494
Protein Weight	32kDa
Shipping	Ice bag