

Anti-Ppap2B Rabbit Polyclonal Antibody For Wb, Elisa - O14495

Item Number: CM0020714

Introduction

Rabbit polyclonal antibody targeting human PPAP2B/PLPP3, suited for WB and ELISA. Reacts with human, mouse, and rat. Ideal for studying lipid phosphate phosphatase in vascular biology and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	PPAP2B Rabbit pAb
Remarks/Alias	LPP3; VCIP; Dri42; PAP2B; PPAP2B
Species	Human
GenelD (Human)	8613
GenelD	8613
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14495
Protein Weight	35kDa
Shipping	Ice bag