


Anti-Ppar α Rabbit Monoclonal Antibody For Wb And Elisa - Q07869

Item Number: CM0011647



Introduction

High-quality PPAR α rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat PPAR α . Ideal for studying lipid metabolism, peroxisome proliferation, and nuclear receptor signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PPAR α Rabbit mAb
Remarks/Alias	PPAR; NR1C1; hPPAR; PPARalpha; PPAR-alpha; PPAR α
Species	Human
GeneID (Human)	5465
GeneID	5465
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07869
Protein Weight	52kDa
Shipping	Ice bag