

Anti-Ppara Rabbit Polyclonal Antibody For Wb - Q07869

Item Number: CM0019605



Introduction

Rabbit polyclonal antibody targeting human PPARA, validated for Western blot. Cross-reacts with mouse. Essential tool for lipid metabolism and nuclear receptor research involving PPAR α .

[Learn More](#)

Parameter	Value
Product Name	PPARA Rabbit pAb
Remarks/Alias	PPAR; NR1C1; hPPAR; PPARalpha; PPARA
Species	Human
Gene ID	5465
Immunogen	Synthetic Peptide corresponding to a sequence within amino acids 201-300 of human PPARA (NP_001001928.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Mouse
SWISS	Q07869
Protein Weight	52kDa
Shipping	Ice bag