


## Anti-Ppard Polyclonal Antibody For Wb, Ihc-P, Elisa - Q03181

Item Number: CM0029820



### Introduction

Rabbit polyclonal PPARD antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets the peroxisome proliferator-activated receptor delta involved in energy metabolism and fatty acid oxidation. ~50kDa.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PPARD Rabbit pAb                                                                                                    |
| Remarks/Alias    | FAAR; NUC1; NUCI; NR1C2; NUCII; PPARB; PPARD                                                                        |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 5467                                                                                                                |
| GeneID           | 5467                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 102-361 of human PPARD (NP_803184.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q03181                                                                                                              |
| Protein Weight   | 50kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |