


Anti-Ppard Rabbit Polyclonal Antibody For Wb And Elisa - Q03181

Item Number: CM0026558



Introduction

Rabbit polyclonal antibody targeting human PPARG, validated for Western blot and ELISA, with cross-reactivity to rat. An essential tool for metabolic and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PPARG Rabbit pAb
Remarks/Alias	FAAR; NUC1; NUCI; NR1C2; NUCI; PPARG
Species	Human
GeneID (Human)	5467
GeneID	5467
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 101-200 of human PPARG (NP_006229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q03181
Protein Weight	50kDa
Shipping	Ice bag