

Anti-Pparγ Rabbit Polyclonal Antibody For IHC-P, IF/ICC, Elisa - P37231

Item Number: CM0013399



Introduction

PPAR γ Rabbit Polyclonal Antibody validated for IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets nuclear receptor PPAR γ , a key regulator of adipogenesis, metabolism, and inflammation. Ideal for metabolic disease and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PPAR γ Rabbit pAb
Remarks/Alias	GLM1; CIMT1; NR1C3; PPARG1; PPARG2; PPARG5; PPARgamma; PPAR γ
Species	Human
GeneID (human)	5468
GeneID	5468
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PPAR γ (NP_056953.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37231
Protein Weight	58kDa
Shipping	Ice bag