

Anti-Ppat Polyclonal Antibody For Wb, If/Icc, Elisa - Q06203

Item Number: CM0028408

Introduction

Rabbit polyclonal anti-PPAT antibody for WB, IF/ICC, ELISA. Recognizes human PPAT (glutamine PRPP amidotransferase), cross-reacts with mouse and rat. ~57 kDa, ideal for purine biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	PPAT Rabbit pAb
Remarks/Alias	GPAT; PRAT; ATASE; PPAT
Species	Human
GeneID (Human)	5471
GeneID	5471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PPAT (NP_002694.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06203
Protein Weight	57kDa
Shipping	Ice bag