


Anti-Ppbb Polyclonal Antibody For Wb, Elisa - P02775

Item Number: CM0022871



Introduction

High-quality anti-PPBP (CXCL7) rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects human platelet basic protein and chemokine variants. Ideal for inflammation, leukocyte activation, and antimicrobial research.

[Learn More](#)

Parameter	Value
Product Name	PPBP Rabbit pAb
Remarks/Alias	PBP; TC1; TC2; TGB; LDGF; MDGF; TGB1; B-TG1; CTAP3; CXCL7; NAP-2; SCYB7; THBGB; LA-PF4; THBGB1; Beta-TG; CTAPIII; CTAP-III; PPBP
Species	Human
GeneID (Human)	5473
GeneID	5473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-128 of human PPBP (NP_002695.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02775
Protein Weight	14kDa
Shipping	Ice bag