


Anti-Ppcs Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Hab8

Item Number: CM0002927



Introduction

Rabbit monoclonal antibody against human PPCS (Phosphopantothenate--cysteine ligase), validated for WB, IHC-P, and ELISA. Essential for studying coenzyme A biosynthesis from vitamin B5.

[Learn More](#)

Parameter	Value
Product Name	PPCS Rabbit mAb
Remarks/Alias	CMD2C; PPCS
Species	Human
GeneID (human)	79717
GeneID	79717
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 212-311 of human PPCS (Q9HAB8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9HAB8
Protein Weight	34kDa
Shipping	Ice bag