

## Anti-Ppfibp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86W92

Item Number: CM0019425

### Introduction

Rabbit polyclonal antibody against human PPFIBP1 (liprin-beta-1) for Western blot and ELISA. Validated in human, mouse, and rat samples. Ideal for cell adhesion and neurological disorder research.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PPFIBP1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | L2; SGT2; hSGT2; hSgt2p; NEDSMBA; PPFIBP1                                                                             |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 8496                                                                                                                  |
| Gene ID          | 8496                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 846-1005 of human PPFIBP1 (NP_003613.3) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q86W92                                                                                                                |
| Protein Weight   | 114kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                               |