


Anti-Pphln1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ney8

Item Number: CM0024552



Introduction

Rabbit polyclonal antibody against human PPHLN1 (Periphilin-1, CR, HSPC206, HSPC232). Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for HUSH complex and epigenetic silencing studies.

[Learn More](#)

Parameter	Value
Product Name	PPHLN1 Rabbit pAb
Remarks/Alias	CR, HSPC206, HSPC232, PPHLN1
Species	Human
Gene ID (Human)	51535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 195-255 of human PPHLN1 (NP_958846.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEY8
Protein Weight	53kDa
Shipping	Ice bag