

Anti-Ppia Rabbit Polyclonal Antibody For Wb, Elisa - P62937

Item Number: CM0016656



Introduction

Rabbit polyclonal antibody targeting human PPIA (Cyclophilin A). Validated for WB and ELISA; cross-reacts with mouse and rat. Useful for studying protein folding, inflammation, apoptosis, and viral infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PPIA Rabbit pAb
Remarks/Alias	CYPA; CYPH; HEL-S-69p; PPIA
Species	Human
GeneID (Human)	5478
GeneID	5478
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-165 of human PPIA (NP_066953.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62937
Protein Weight	18kDa
Shipping	Ice bag