

Anti-Ppic Rabbit Polyclonal Antibody For Wb, Elisa - P45877

Item Number: CM0027212



Introduction

CamelBio's PPIC rabbit polyclonal antibody, specific to human cyclophilin C, enables Western blot and ELISA detection. Designed for researchers investigating protein folding and isomerase activity.

[Learn More](#)

Parameter	Value
Product Name	PPIC Rabbit pAb
Remarks/Alias	CYPC; PPIC
Species	Human
Gene ID (Human)	5480
Gene ID	5480
Immunogen	Recombinant fusion protein containing amino acids 23-212 of human PPIC (NP_000934.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P45877
Protein Weight	23kDa
Shipping	Ice bag