

Anti-Ppie Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Unp9

Item Number: CM0026161



Introduction

High-quality Anti-PPIE rabbit polyclonal antibody validated for WB, IP, ELISA. Reacts with human and mouse. Peptidyl-prolyl cis-trans isomerase E involved in pre-mRNA splicing and proline isomerization.

[Learn More](#)

Parameter	Value
Product Name	PPIE Rabbit pAb
Remarks/Alias	CYP33, CYP-33, PPIE
Species	Human
GeneID (human)	10450
GeneID	10450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human PPIE (NP_006103.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNP9
Protein Weight	33kDa
Shipping	Ice bag