


Anti-Ppif (Cyclophilin D) Polyclonal Antibody For Wb And Elisa - P30405

Item Number: CM0027246



Introduction

Rabbit polyclonal antibody against human PPIF (Cyclophilin D), validated for Western blotting (WB) and ELISA. Cross-reacts with mouse. Suitable for mitochondrial permeability transition pore and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PPIF Rabbit pAb
Remarks/Alias	CYP3, CypD, CyP-M, Cyp-D
Species	Human
GeneID (Human)	10105
GeneID	10105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of human PPIF (NP_005720.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P30405
Protein Weight	22kDa
Shipping	Ice bag