


Anti-Ppif Rabbit Polyclonal Antibody For Wb And Elisa - P30405

Item Number: CM0030568



Introduction

Rabbit polyclonal antibody against human PPIF (Cyclophilin D), suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Targets mitochondrial matrix protein involved in permeability transition pore regulation and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	PPIF Rabbit pAb
Remarks/Alias	CYP3; CypD; CyP-M; Cyp-D; PPIF
Species	Human
GeneID (human)	10105
GeneID	10105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-207 of human PPIF (NP_005720.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30405
Protein Weight	22kDa
Shipping	Ice bag