

Anti-Ppl Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60437

Item Number: CM0030154



Introduction

Rabbit polyclonal antibody targeting periplakin (PPL), validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Suitable for keratinocyte cornified envelope and desmosome research.

[Learn More](#)

Parameter	Value
Product Name	PPL Rabbit pAb
Remarks/Alias	PPL
Species	Human
GeneID (Human)	5493
GeneID	5493
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60437
Protein Weight	205kDa
Shipping	Ice bag