

Anti-Ppm1A Monoclonal Antibody For Wb, Elisa - P35813

Item Number: CM0002212



Introduction

Anti-PPM1A Rabbit Monoclonal Antibody validated for WB and ELISA targets human PPM1A/PP2C-alpha, a phosphatase in TGF-beta and NF-kB pathways. Suitable for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PPM1A Rabbit mAb
Remarks/Alias	PP2CA; PP2Calpha; PP2C-ALPHA; PPM1A
Species	Human
GeneID (Human)	5494
GeneID	5494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-382 of human PPM1A (NP_066283.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35813
Protein Weight	42kDa
Shipping	Ice bag