

Anti-Ppm1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35813

Item Number: CM0025879



Introduction

Rabbit polyclonal anti-PPM1A antibody for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat PPM1A (PP2C α , 42 kDa). Ideal for TGF-beta and NF-kB signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PPM1A Rabbit pAb
Remarks/Alias	PP2CA; PP2Calpha; PP2C-ALPHA; PPM1A
Species	Human
GeneID (Human)	5494
GeneID	5494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 161-370 of human PPM1A (NP_005029.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35813
Protein Weight	42kDa
Shipping	Ice bag