

Anti-Ppm1F Rabbit Polyclonal Antibody For Wb, Elisa - P49593

Item Number: CM0015897

Introduction

High specificity rabbit polyclonal antibody targeting human PPM1F (CaM-kinase phosphatase), suitable for Western blot and ELISA. Ideal for apoptosis and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PPM1F Rabbit pAb
Remarks/Alias	CAMKP; FEM-2; POPX2; hFEM-2; CaMKPase; PPM1F
Species	Human
GeneID (Human)	9647
GeneID	9647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human PPM1F (NP_055449.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49593
Protein Weight	50kDa
Shipping	Ice bag