

Anti-Ppm1J Rabbit Polyclonal Antibody For Wb And Elisa - Q5Jr12

Item Number: CM0018296

Introduction

Rabbit polyclonal antibody against human PPM1J (PP2C-zeta) for WB and ELISA. Validated with recombinant immunogen. Cross-reacts with mouse. UniProt Q5JR12. Ideal for research on protein phosphatase 1J, a Ser/Thr phosphatase involved in cellular signaling. Suitable for Western blot and ELISA.

[Learn More](#)

Parameter	Value
Product Name	PPM1J Rabbit pAb
Remarks/Alias	PP2CZ, PPP2CZ, PP2Czeta, PP2C-zeta, PPM1J
Species	Human
GeneID	333926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-260 of human PPM1J (NP_005158.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5JR12
Protein Weight	55kDa
Shipping	Ice bag