


Anti-Ppme1 Polyclonal Antibody For Wb And Elisa - Q9Y570

Item Number: CM0017148



Introduction

Rabbit polyclonal antibody against PPME1 (PME-1), suitable for WB and ELISA. Reacts with human, mouse, and rat. Detects protein phosphatase methylesterase 1, a key regulator of PP2A methylation and activity.

[Learn More](#)

Parameter	Value
Product Name	PPME1 Rabbit pAb
Remarks/Alias	PME-1; ABDH19; PPME1
Species	Human
GeneID (Human)	51400
GeneID	51400
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human PPME1 (NP_057231.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y570
Protein Weight	42kDa
Shipping	Ice bag