

Anti-Ppp1Ca Polyclonal Antibody For Wb, If/Icc, Elisa - P62136

Item Number: CM0018989



Introduction

Rabbit polyclonal antibody targeting human PPP1CA (PP1alpha catalytic subunit). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. MW 38kDa. Ideal for phosphatase signaling and MAPK pathway research.

[Learn More](#)

Parameter	Specification
Product Name	PPP1CA Rabbit pAb
Remarks/Alias	PP1A, PP-1A, PPP1A, PP1alpha, PPP1CA
Species	Human
Gene ID	5499 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human PPP1CA (NP_002699.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62136
Protein Weight	38kDa
Shipping	Ice bag