

Anti-Ppp1Ca Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62136

Item Number: CM0006438



Introduction

High-quality rabbit monoclonal antibody targeting human PPP1CA (PP1-alpha catalytic subunit). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for phosphatase research in cell cycle, transcription, and signaling.

[Learn More](#)

Parameter	Value
Product Name	PPP1CA Rabbit mAb
Remarks/Alias	PP1A; PP-1A; PPP1A; PP1alpha; PPP1CA
Species	Human
GeneID (Human)	5499
GeneID	5499
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 230-330 of human PPP1CA (NP_002699.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62136
Protein Weight	38kDa
Shipping	Ice bag