

Anti-Ppp1Ca Rabbit Polyclonal Antibody For Wb - P62136

Item Number: CM0016319



Introduction

Rabbit polyclonal antibody targeting human PPP1CA (PP1-alpha), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. For studies on phosphatase signaling, cell cycle, and metabolism.

[Learn More](#)

Parameter	Value
Product Name	PPP1CA Rabbit pAb
Remarks/Alias	PP1A; PP-1A; PPP1A; PP1alpha; PPP1CA
Species	Human
Gene ID (Human)	5499
Gene ID	5499
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-330 of human PPP1CA (NP_002699.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62136
Protein Weight	38kDa
Shipping	Ice bag