


Anti-Ppp1R12A Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O14974

Item Number: CM0014276



Introduction

Rabbit polyclonal antibody targeting PPP1R12A (MYPT1), validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Key regulator of myosin phosphatase and PPP1C complex.

[Learn More](#)

Parameter	Value
Product Name	PPP1R12A Rabbit pAb
Remarks/Alias	MBS; GUBS; M130; MYPT1; PPP1R12A
Species	Human
GeneID	4659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PPP1R12A (NP_002471.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14974
Protein Weight	115kDa
Shipping	Ice bag