

Anti-Ppp1R12B Rabbit Polyclonal Antibody For Wb, Elisa - O60237

Item Number: CM0019990

Introduction

Rabbit polyclonal antibody targeting PPP1R12B (MYPT2) validated for WB and ELISA. Reacts with human; ideal for studying myosin phosphatase regulation and cytoskeletal organization.

[Learn More](#)

Parameter	Value
Product Name	PPP1R12B Rabbit pAb
Remarks/Alias	MYPT2; PP1bp55; PPP1R12B
Species	Human
GeneID (human)	4660
GeneID	4660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 660-930 of human PPP1R12B (NP_002472.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60237
Protein Weight	110kDa
Shipping	Ice bag