

Anti-Ppp1R13L Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wuf5

Item Number: CM0014614



Introduction

Rabbit polyclonal antibody targeting human PPP1R13L (iASPP, RelA-associated inhibitor), suitable for WB and ELISA. Uses include studying apoptosis, NF-kappa-B-dependent inflammatory response, and p53 transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	PPP1R13L Rabbit pAb
Remarks/Alias	RAI; RAI4; IASPP; NKIP1; PPP1R13L
Species	Human
GeneID (Human)	10848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 589-828 of human PPP1R13L (NP_006654.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WUF5
Protein Weight	89kDa
Shipping	Ice bag