

Anti-Ppp1R14A (Cpi-17) Rabbit Polyclonal Antibody For Wb And Elisa - Q96A00

Item Number: CM0016816

Introduction

Rabbit polyclonal antibody targeting PPP1R14A/CPI-17, a phosphorylation-dependent inhibitor of PP1 implicated in smooth muscle regulation. Validated for Western blot and ELISA, cross-reactive with human, mouse, and rat. Ideal for cellular signaling and contraction research.

[Learn More](#)

Parameter	Value
Product Name	PPP1R14A Rabbit pAb
Remarks/Alias	CPI17; CPI-17; PPP1INL; PPP1R14A
Species	Human
Gene ID (Human)	94274
Gene ID	94274
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human PPP1R14A (NP_150281.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96A00
Protein Weight	17kDa
Shipping	Ice bag