

## Anti-Ppp1R14B Rabbit Polyclonal Antibody For Wb, Elisa - Q96C90

Item Number: CM0016839

### Introduction

Rabbit polyclonal antibody against human PPP1R14B. Validated for WB and ELISA, reacts with human, mouse, and rat. Phosphorylation-sensitive PP1 inhibitor; ~16 kDa. Ideal for phosphatase signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PPP1R14B Rabbit pAb                                                                                                  |
| Remarks/Alias    | PNG; PHI-1; PLCB3N; SOM172; PPP1R14B                                                                                 |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 26472                                                                                                                |
| GeneID           | 26472                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 68-147 of human PPP1R14B (NP_619634.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | Q96C90                                                                                                               |
| Protein Weight   | 16kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |