

Anti-Ppp1R15A Rabbit Polyclonal Antibody (Gadd34) For Wb, If/Icc, Elisa - O75807

Item Number: CM0022129



Introduction

Anti-PPP1R15A rabbit polyclonal antibody (GADD34) for WB, IF/ICC, ELISA. Validated in human, mouse, rat. Targets stress-inducible phosphatase regulatory subunit involved in translation recovery and autophagy.

[Learn More](#)

Parameter	Value
Product Name	PPP1R15A Rabbit pAb
Remarks/Alias	GADD34; PPP1R15A
Species	Human
GeneID (Human)	23645
GeneID	23645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 375-674 of human PPP1R15A (NP_055145.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75807
Protein Weight	73kDa
Shipping	Ice bag