

Anti-Ppp1R3B Rabbit Polyclonal Antibody For Wb, Elisa - Q86Xi6

Item Number: CM0018198



Introduction

Rabbit polyclonal antibody targeting human PPP1R3B (glycogen-targeting subunit GL). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Suitable for glycogen metabolism and insulin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PPP1R3B Rabbit pAb
Remarks/Alias	GL; PTG; PPP1R4; PPP1R3B
Species	Human
Gene ID (Human)	79660
Gene ID	79660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human PPP1R3B (NP_078883.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86Xi6
Protein Weight	33kDa
Shipping	Ice bag