

Anti-Ppp1R3C Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9Uqk1

Item Number: CM0021964



Introduction

High-quality rabbit polyclonal antibody against human PPP1R3C (PTG), validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for glycogen metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PPP1R3C Rabbit pAb
Remarks/Alias	PTG, PPP1R5, PPP1R3C
Species	Human
GeneID (Human)	5507
GeneID	5507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human PPP1R3C (NP_005389.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQK1
Protein Weight	36kDa
Shipping	Ice bag