


Anti-Ppp1R7 Rabbit Polyclonal Antibody For Wb, Elisa - Q15435

Item Number: CM0016111



Introduction

This PPP1R7 rabbit polyclonal antibody detects endogenous levels of protein phosphatase 1 regulatory subunit 7 in human and mouse samples. Validated for Western blot (WB) and ELISA, it targets the 281-360 amino acid region. Ideal for IVD research use.

[Learn More](#)

Parameter	Value
Product Name	PPP1R7 Rabbit pAb
Remarks/Alias	SDS22; PPP1R7
Species	Human
Gene ID	5510
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-360 of human PPP1R7 (NP_002703.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15435
Protein Weight	42kDa
Shipping	Ice bag