

## Anti-Ppp1R9B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Sb3

Item Number: CM0023436

### Introduction

Rabbit polyclonal antibody against human PPP1R9B (Neurabin-2/Spinophilin), cross-reactive with human and mouse. Suitable for WB, IF/ICC, ELISA. Ideal for studying synaptic transmission, dendritic spines, and actin cytoskeleton interactions.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PPP1R9B Rabbit pAb                                                                                                           |
| Remarks/Alias    | Spn; SPINO; PPP1R6; PPP1R9; PPP1R9B                                                                                          |
| Species          | Human                                                                                                                        |
| Gene ID (Human)  | 84687                                                                                                                        |
| Gene ID          | 84687                                                                                                                        |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PPP1R9B (NP_115984.3). |
| Source           | Rabbit                                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                                        |
| Application      | WB, IF/ICC, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                                 |
| SWISS            | Q96SB3                                                                                                                       |
| Protein Weight   | 89kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |