

Anti-Ppp2R1A/Ppp2R2A/Ppp2R1B Rabbit Polyclonal Antibody For Wb, Elisa - P30153 / P63151 / P30154

Item Number: CM0022302



Introduction

High-quality rabbit polyclonal antibody targeting PPP2R1A, PPP2R2A, and PPP2R1B. Suitable for Western blotting and ELISA in human, mouse, and rat. Investigate PP2A complex function in signaling, mitosis, and transcription.

[Learn More](#)

Parameter	Value
Product Name	PPP2R1A/PPP2R2A/PPP2R1B Rabbit pAb
Remarks/Alias	PPP2R1A/PPP2R2A/PPP2R1B
Species	Human
GeneID (Human)	5518/5520/5519
GeneID	5518/5520/5519
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human PPP2R1A/PPP2R2A/PPP2R1BPPP2R1A (NP_055040.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30153/P63151/P30154
Shipping	Ice bag