

Anti-Ppp2R2B/Ppp2R2C Polyclonal Antibody For Wb, Elisa - Q00005 / Q9Y2T4

Item Number: CM0021873

Introduction

Rabbit polyclonal antibody targeting PPP2R2B/PPP2R2C, the regulatory subunits of protein phosphatase 2A. Validated for WB and ELISA, reacts with human, mouse, and rat samples. Ideal for cell signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PPP2R2B/PPP2R2C Rabbit pAb
Remarks/Alias	PPP2R2B/PPP2R2C
Species	Human
Gene ID (Human)	5521/5522
Gene ID	5521/5522
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 108-220 of human PPP2R2B/PPP2R2C (NP_858060.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00005 / Q9Y2T4
Shipping	Ice bag