

Anti-Ppp2R3B Polyclonal Antibody For Wb And Elisa - Q9Y5P8

Item Number: CM0028716



Introduction

Rabbit polyclonal antibody targeting human PPP2R3B, a regulatory subunit of PP2A. Validated for WB and ELISA, with cross-reactivity to human and mouse. Suitable for studying PP2A function and nuclear localization.

[Learn More](#)

Parameter	Value
Product Name	PPP2R3B Rabbit pAb
Remarks/Alias	PR48; PR70; NYREN8; PPP2R3L; PPP2R3LY; PPP2R3B
Species	Human
GeneID (Human)	28227
GeneID	28227
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 356-575 of human PPP2R3B (NP_037371.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5P8
Protein Weight	65kDa
Shipping	Ice bag