

Anti-Ppp3Cc Polyclonal Antibody For Wb, Elisa - P48454

Item Number: CM0027973



Introduction

Rabbit polyclonal antibody targeting PPP3CC (calcineurin A gamma), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for calcium signaling and testis-specific phosphatase research.

[Learn More](#)

Parameter	Value
Product Name	PPP3CC Rabbit pAb
Remarks/Alias	CNA3; CALNA3; PP2Bgamma; PPP3CC
Species	Human
Gene ID	5533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 373-512 of human PPP3CC (NP_005596.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48454
Protein Weight	58kDa
Shipping	Ice bag