

Anti-Ppp3R2 Polyclonal Antibody For Wb, Elisa - Q96Lz3

Item Number: CM0028562

Introduction

PPP3R2 rabbit polyclonal antibody for WB and ELISA, specific to human, mouse, and rat. Recognizes calcineurin B type 2, a testis-specific regulatory subunit involved in calcium-dependent phosphatase activity.

[Learn More](#)

Parameter	Value
Product Name	PPP3R2 Rabbit pAb
Remarks/Alias	PPP3RL; PPP3R2
Species	Human
GeneID (Human)	5535
GeneID	5535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human PPP3R2 (NP_671709.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96LZ3
Protein Weight	20kDa
Shipping	Ice bag