

Anti-Ppp4C Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P60510

Item Number: CM0013278

Introduction

Rabbit polyclonal antibody against human PPP4C, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for DNA repair, apoptosis, and signaling research. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	PPP4C Rabbit pAb
Remarks/Alias	PP4; PPX; PP-X; PP4C; PPH3; PPP4; PPP4C
Species	Human
GeneID (Human)	5531
GeneID	5531
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-307 of human PPP4C (NP_002711.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60510
Protein Weight	35kDa
Shipping	Ice bag