

## Anti-Ppp4R1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q8Tf05

Item Number: CM0030243

### Introduction

Anti-PPP4R1 rabbit polyclonal antibody validated for WB, IP, ELISA in human and mouse samples. Targets serine/threonine-protein phosphatase 4 regulatory subunit 1, a key regulator of NF-kappa-B signaling and cell division. Ideal for IVD raw material sourcing.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PPP4R1 Rabbit pAb                                                                                                       |
| Remarks/Alias    | MEG1; PP4R1; PP4(Rmeg); PPP4R1                                                                                          |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 9989                                                                                                                    |
| GeneID           | 9989                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 324-647 of human PPP4R1 (NP_001035847.1). |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IP, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                            |
| SWISS            | Q8TF05                                                                                                                  |
| Protein Weight   | 107kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |