

Anti-Ppp5C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P53041

Item Number: CM0017136



Introduction

High-quality rabbit polyclonal antibody targeting human PPP5C (PP5/PPT), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat; ideal for studying DNA damage response, apoptosis, and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	PPP5C Rabbit pAb
Remarks/Alias	PP5, PPT, PPP5, PPP5C
Species	Human
GeneID (Human)	5536
GeneID	5536
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PPP5C (NP_006238.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53041
Protein Weight	57kDa
Shipping	Ice bag