

Anti-Ppp6C Rabbit Monoclonal Antibody For Wb, Elisa - 000743

Item Number: CM0009070

Introduction

Rabbit monoclonal antibody against PPP6C (Serine/threonine-protein phosphatase 6 catalytic subunit). Suitable for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying cell cycle regulation and innate immunity pathways.

[Learn More](#)

Parameter	Value
Product Name	PPP6C Rabbit mAb
Remarks/Alias	PP6; PP6C; PPP6C
Species	Human
GeneID	5537
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00743
Protein Weight	35kDa
Shipping	Ice bag