


Anti-Ppt1 Rabbit Polyclonal Antibody For Wb, Elisa - P50897

Item Number: CM0018181



Introduction

Rabbit polyclonal antibody to human PPT1, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for studying palmitoyl-protein thioesterase function in lysosomal storage diseases.

[Learn More](#)

Parameter	Value
Product Name	PPT1 Rabbit pAb
Remarks/Alias	PPT; CLN1; INCL; PPT1
Species	Human
GeneID (Human)	5538
GeneID	5538
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-306 of human PPT1 (NP_000301.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50897
Protein Weight	34kDa
Shipping	Ice bag