

Anti-Pqbp1 Rabbit Monoclonal Antibody For Wb, Elisa - O60828

Item Number: CM0002263



Introduction

Rabbit monoclonal antibody against human PQBP1 (30 kDa), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying pre-mRNA splicing, transcription, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	PQBP1 Rabbit mAb
Remarks/Alias	SHS; MRX2; MRX55; MRX53; MRX58; NPW38; RENS1; PQBP1
Species	Human
GeneID	10084
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human PQBP1 (NP_001027553.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60828
Protein Weight	30kDa
Shipping	Ice bag