

Anti-Pqbp1 Rabbit Polyclonal Antibody For Wb And Elisa - O60828

Item Number: CM0027563



Introduction

Rabbit polyclonal antibody against human PQBP1, validated for Western blotting (WB) and ELISA. Reacts with human, mouse, and rat. PQBP1 is an intrinsically disordered scaffold protein involved in pre-mRNA splicing, transcription regulation, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	PQBP1 Rabbit pAb
Remarks/Alias	SHS; MRX2; MRX55; MRX53; MRX58; NPW38; RENS1; PQBP1
Species	Human
GeneID	10084 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human PQBP1 (NP_005701.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60828
Protein Weight	30kDa
Shipping	Ice bag