

Anti-Praf2 Rabbit Polyclonal Antibody For Wb, Elisa - O60831

Item Number: CM0021037

Introduction

Rabbit polyclonal antibody targeting human PRAF2 (JM4, Yip6a), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis and vesicular traffic studies.

[Learn More](#)

Parameter	Value
Product Name	PRAF2 Rabbit pAb
Remarks/Alias	JM4; Yip6a; PRAF2
Species	Human
GeneID (Human)	11230
GeneID	11230
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-50 of human PRAF2 (NP_009144.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60831
Protein Weight	19kDa
Shipping	Ice bag